


AWARD-WINNING & INNOVATIVE

Oxford SpheriChrome

One product = multiple shades



Dye & Pigment free nano-hybrid composite with spherical fillers produces a chameleon effect and adapts to match the patient's specific dentin shade. Flowable version, Oxford SpheriChrome Flow, available with an optimal thixotropic consistency.



The shape and distribution of **Oxford SpheriChrome's** filler particles allows white light to pass through the composite and combines with the colors of the patients surrounding dentin/enamel to reflect a patient specific shade. In cases with severe tooth discoloration or extended Class III and IV cavities; **Oxford SpheriBlock** is available for masking and blocking.

CLINICAL CASES

Performed by Dr. Ross Nash, DDS & Dr. Tyler Wurmlinger, DDS

Case 1:

Patient with wear and occlusal carries in mandibular left first and second molars. Restoration performed with Oxford SpheriChrome Flow.

Case 2:

Patient asked to have the central diastema closed for a more "normal-looking" smile. All four incisors were directly veneered with Oxford SpheriChrome. Oxford SpheriBlock was needed at the interproximal space between the central incisors.



Photos by Dr. Ross Nash




Oxford SpheriChrome produces excellent restorations with less effort and fits easily into the everyday office procedure while delivering on claims to virtually eliminate shade matching. It saves time, yields improved restorations, reduces inventory, and increases patient flow in the office.

Oxford SpheriChrome's Superior Physical Properties

Tooth shade isn't the only important factor to consider; strength, surface finish, and shrinkage impacts the overall life and performance of the restoration.

- ✓ Smoother Surface roughness = better aesthetics & wear; higher gloss less discoloration
- ✓ Higher Strength = higher wear resistance to withstand chewing forces
- ✓ Lower shrinkage = better marginal quality & reduced chance for micro leakage.



Ordering Information

Oxford SpheriChrome

- 10-007 Oxford SpheriChrome 4 gr. Syringe
- 10-011 Oxford SpheriChrome 20 x 0.25g capsules
- 11-004 Oxford SpheriChrome Flow 2 x 1 ml Syringes + 8 needle tips

Oxford SpheriBlock

- 10-010 Oxford SpheriBlock 4 gr Syringe
- 10-012 Oxford SpheriBlock 20 x 0.25g capsules
- 11-005 SpheriBlock Flow 2 x 1 ml Syringes + 8 needle tips

5 Key Reasons Oxford SpheriChrome is a superior material

1. Hyper-Nano technology
2. Spherical fillers
3. Pigment free
4. Light optical effect for a great match with natural tooth substance
5. Low shrinkage



Oxford Scientific
581 Inner Circle, The Villages, FL 32162
Orders@OxfordScientificNA.com
Ph: 908.553.3196

www.OxfordScientificNA.com

Oxford Distributor Partner: